

## Care Instructions - GLING™ muscle roller

The GLING™ muscle roller is made to blend in to a living room type environment as much as serving as a pre-workout tool or a part of daily recovery routine after either a tough workout or a hard day at work. Crafted from one single, hand-stitched piece of leather, this products has solid wood ends with stainless steel inserted end-caps that results in a classic design that can be appreciated both during use and while on display.

To ensure the longevity of this handcrafted product, we recommend occasional maintenance. Approximately every 6-12 months, apply a leather conditioner, preferably one made from beeswax, to keep the leather supple and looking its best. Additionally, avoid storing the muscle roller on sharp edges or objects to prevent unsightly dents and maintain the pristine condition of the leather over time.

